

# The DISTAL GLIDE™ - effective distalization



- Pre-manufactured distal glide spring for consistent precision
- CNA Beta III/TM Nickel-free combines wires which offers excellent resiliency when compared to Stainless Steel
- Easy oral hygiene
- Easily removed for activation outside the patients mouth
- Full control over torque, tip and rotation
- Clear acrylic Nance-plate for better visual control of hygiene and possible infections

The DISTAL GLIDE™ (modified Pendulum) has proved to be an incredibly efficient and predictable source of non-compliant Class II correction. Our aim was to develop the simplest, easy to handle appliance that would allow the patient to maintain good oral hygiene, improve patient comfort, ease of appliance placement and activation, and obtain the most effective distalization.



**Quality** - adenta engineered an appliance to provide increased stability and improved overall response. The distalizing springs are of high quality Nickel-free Beta Titanium.

**Pre-engineered** - The distalizing springs are fabricated by machine with absolute precision. Each spring is identical to each other, and provides the practitioner with consistent forces.

**Removable** - To reactivate the old Pendulum appliances intraorally was cumbersome and was uncomfortable for the patient. The Distal Glide springs are inserted into tubes located inside the acrylic plate and into the lingual sheath, this provides a sturdy support, and can easily be removed for reactivation.

**Hygienic** - As the distalizing springs are easily removed, this provides easier access for the hygienist.

## Distal Glide™ order info

Item	Article Description
DP-DGP-Info	ADENTA Distal Glide Information
DP-DGP-PR	Distal Glide - Office Kit, 10 distalizing
DP-DGP-ST	Distal Glide - Starter Kit, 5 distalizing