

Tooth	Torque	Tip	Item # .018	Item # .022
11 - Central	12°	5°	105-11	155-11
12 - Lateral	8°	9°	105-12	155-12
13 - Cuspid	-2°	9°	105-13	155-13
13 - Cuspid w hook	-2°	9°	105-13/H	155-13/H
14 - 1. Bicuspid	-7°	0°	105-14/15	155-14/15
14 - 1. Bicuspid w hook	-7°	0°	105-14/15/H	155-14/15/H
15 - 2. Bicuspid	-7°	0°	105-14/15	155-14/15
15 - 2. Bicuspid w hook	-7°	0°	105-14/15/H	155-14/15/H

Item # .022	Item # .018	Tip	Torque	Tooth
155-21	105-21	5°	12°	Central - 21
155-22	105-22	9°	8°	Lateral - 22
155-23	105-23	9°	-2°	Cuspid - 23
155-23/H	105-23/H	9°	-2°	Cuspid w hook - 23
155-24/25	105-24/25	0°	-7°	1. Bicuspid - 24
155-24/25/H	105-24/25/H	0°	-7°	1. Bicuspid w hook - 24
155-24/25	105-24/25	0°	-7°	2. Bicuspid - 25
155-24/25/H	105-24/25/H	0°	-7°	2. Bicuspid w hook - 25

Tooth	Torque	Tip	Item # .018	Item # .022
41 - Anterior	-1°	0°	105-41/42	155-41/42
42 - Anterior	-1°	0°	105-41/42	155-41/42
43 - Cuspid	-11°	7°	105-43	155-43
43 - Cuspid w hook	-11°	7°	105-43/H	155-43/H
44 - 1. Bicuspid	-17°	0°	105-44	155-44
44 - 1. Bicuspid w hook	-17°	0°	105-44/H	155-44/H
45 - 2. Bicuspid	-22°	0°	105-45	155-45
45 - 2. Bicuspid w hook	-22°	0°	105-45/H	155-45/H

Item # .022	Item # .018	Tip	Torque	Tooth
155-31/32	105-31/32	0°	-1°	Anterior - 31
155-31/32	105-31/32	0°	-1°	Anterior - 32
155-33	105-33	7°	-11°	Cuspid - 33
155-33/H	105-33/H	7°	-11°	Cuspid w hook - 33
155-34	105-34	0°	-17°	1. Bicuspid - 34
155-34/H	105-34/H	0°	-17°	1. Bicuspid w hook - 34
155-35	105-35	0°	-22°	2. Bicuspid - 35
155-35/H	105-35/H	0°	-22°	2. Bicuspid w hook - 35